

Recombinant Human Junctional Adhesion Molecule B, JAM-B, CD322 (C-Fc)



Catalog No: #AP75554

Orders: order@signalwayantibody.com

Package Size: #AP75554-1 10ug #AP75554-2 50ug #AP75554-3 500ug #AP75554-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Junctional Adhesion Molecule B, JAM-B, CD322 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P57087
Uniprot	P57087
GeneID	58494;
Calculated MW	50.3kDa
Target Sequence	FSAPKDQQVVTAVEYQEAILACKTPKKTVSSRLEWKKLGRSVSFVYYQQLQGDFKNRAEMIDFNIRIKNVTRS DAGKYRCEVSAPSEQGQNLEEDTVTLEVLVAPAVPSCVPEVPSSALSGTVVELRCQDKEGNPAPEYTWFKDGIR LLENPRLGSQSTNSSYTMNTKTGTLQFNVTVSKLDTGEYSCEARNVSGYRRCPGKRMQVDDLNVDDIEGRMD EPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNA KTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTK NQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLTVDKSRWQQGNVFCFSVMHEAL HNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4 .
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only