

Recombinant Human ADAM-like Protein Decysin-1,ADAMDEC1 (C-6His)



Catalog No: #AP75556

Orders: order@signalwayantibody.com

Package Size: #AP75556-1 10ug #AP75556-2 50ug #AP75556-3 500ug #AP75556-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ADAM-like Protein Decysin-1,ADAMDEC1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O15204
Uniprot	O15204
GeneID	27299;
Calculated MW	30.3kDa
Target Sequence	KSPEKEDFLRAQKYIDLYLVLDNAFYKNYNENLTLIRSFVFDVMNLLNVIYNTIDVQVALVGMEIWSGDGKIKVV PSASTTFDNFLRWHSSNLGKKIHDHAQLLSGISFNRRRVGLAASNLCSPSSVAVIEAKKKNNVALVGVMSHE LGHVLGMPDVPFNTKCPGSCVMNQYLSSKFPKDFSTSCRAHFERYLLSQPKCLLQAPIPTNIMTTPVCGNH LLEVGEDCDCGSPKECTNLCCEALTCKLPGTDCGGDAPNHTTEVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mM NaCl,1mM ZnCl2,pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only