

# Recombinant Human Peptidyl-Prolyl Cis-Trans Isomerase FKBP2,FKBP22,FKBP13 (C-6His)



Catalog No: #AP75557

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75557-1 10ug #AP75557-2 50ug #AP75557-3 500ug #AP75557-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Peptidyl-Prolyl Cis-Trans Isomerase FKBP2,FKBP22,FKBP13 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P26885
Uniprot	P26885
GeneID	2286;
Calculated MW	14.3kDa
Target Sequence	ATGAEGKRKRLQIGVKKRVDHCPKSRKGDVLHMHYTGKLEDGTEFDSSLPQNQPFFVSLGTGQVIKGDQGL LGMCEGEKRKLVIPSELGYGERGAPPKIPGGATLVFEVELLKIERRTELVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mm NaCl,pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only