

# Recombinant Human Carboxypeptidase E,CPE (C-6His)

Catalog No: #AP75558

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75558-1 10ug #AP75558-2 50ug #AP75558-3 500ug #AP75558-4 1mg

## Description

Product Name	Recombinant Human Carboxypeptidase E,CPE (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P16870
Uniprot	P16870
GeneID	1363;
Calculated MW	47.1kDa
Target Sequence	RLQQEDGISFEYHRYPELREALVSVWLQCTAISRIYTVGRSFEGRELLVIELSDNPGVHEPGEPEFKYIGNMHG NEAVGRELLIFLAQYLCNEYQKGNETIVNLIHSTRIHIMPSLNPDGFEKAASQPGEKDWDFVGRSNAQQIDLNR NFPDLDRIVYVNEKEGGPNNHLLKNMKKIVDQNTKLAPETKAVIHWIMDIPFVLSANLHGGDLVANYPYDETRS GSAHEYSSPDDAIFQSLARAYSSFNPAMSDPNRPPCRKNDDSSFVDGTTNGGAWYSVPGGMQDFNYLSS NCFEITVELSCEKFPPEETLKTYWEDNKNSLISYLEQIHRGVKGFVRDLQGNPIANATISVEGIDHDVTSKDG YWRLLIPGNYKLTASAPGYLAITKKVAVPYSPAAGVDFELESFSDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mM NaCl,1mM ZnCl <sub>2</sub> ,10%Glycerol,pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only