

Recombinant Human Carboxypeptidase M,CPM (C-6His)

Catalog No: #AP75559

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75559-1 10ug #AP75559-2 50ug #AP75559-3 500ug #AP75559-4 1mg

Description

Product Name	Recombinant Human Carboxypeptidase M,CPM (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P14384
Uniprot	P14384
GeneID	1368;
Calculated MW	47.3kDa
Target Sequence	<p> LDFNYHRQEGMEAFKLTVAQNYSSVTHLSIGKSVKGRNLWVLVVGFRFPKEHRIGIPEFKYVANMHGDETVG RELLHLIDYLVTSDGKDPEITNLINSTRIHIMPSMNPDGFEAVKKPDCYYSIGREYNQYDLNRNFPDAFEYNN VSRQPETVAVMKWLKTETFLSANLHGGALVASYPFDNGVQATGALYSRSLTPDDDFQYLAHTYASRNPNM KKGDECKNKMNFPNGVTNGYSWYPLQGGMQDYNIWAQCFEITLELSCCKYPREEKLPSEFWNNKASLIEYI KQVHLGVKGVFDQNGNPLPNVIVEVQDRKHICPYRTNKYGEYLLLLLPGSYIINVTVPGHDPHITKVIPEKSQ NFSALKKDILLPFQQLDSIPVSNPSCPMIPLYRNLDPDHVDHHHHHH </p>
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mM NaCl,1mM ZnCl ₂ ,pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only