

Recombinant Human Ketohexokinase, KHK (C-6His)

Catalog No: #AP75561



Package Size: #AP75561-1 10ug #AP75561-2 50ug #AP75561-3 500ug #AP75561-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Ketohexokinase, KHK (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P50053
Uniprot	P50053
GeneID	3795;
Calculated MW	33.7kDa
Target Sequence	MEEKQILCVGLVVLVDVISLVDKYPKEDSEIRCLSQRWQRGGNASNSCTILSLLGAPCAFMGSMAPGHVADFVL DDLRRYSVDLRYTVFQTTGSVPIATVIINEASGSRTILYYDRSLPDVSATDFEKVDLTQFKWIHIEGRNASEQVK MLQRIDAHNTRQPPEQKIRVSVEVEKPREELFQLFGYGDVVFVSKDVAKHLGFQSAEEALRGLYGRVRKGAV LVCAWAEEGADALGPDGKLLHSDAFPPPRVVDTLGAGDTFNASVIFLSQGRSVQEALRFGCQVAGKKCGLQ GFDGIVVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 50nM KCl, 10% Glycerol, pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only