

# Recombinant Human Pregnancy-Specific β-1-Glycoprotein 5,PSG5 (C-6His)

Catalog No: #AP75562

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75562-1 10ug #AP75562-2 50ug #AP75562-3 500ug #AP75562-4 1mg

## Description

Product Name	Recombinant Human Pregnancy-Specific β-1-Glycoprotein 5,PSG5 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q15238
Uniprot	Q15238
GeneID	5673;
Calculated MW	35kDa
Target Sequence	QVTIEALPPKVSEGKDVLLLVHNLQPQLAGYIWYKGQLMDLYHYITSYVVDGQINIYGPAYTGRETVYSNASLLI QNVTRDAGSYTLHIKRGDRTRGVTFYFTFNLYLKLKPKPYITINNSKPRENKDVLAFTCPEKSENYTIWWLN GQSLPVS PRVKQPIENRILILPSVTRNETGPYECEIRDRDGGMHSDPVTNLNVLGPDLPISIYPSFTYYRSGENL YLSCFAESNPPAEYFWTINGKFQQSGQKLSIPQITTKHRGLYTCSVRNSATGKESSKSMTEVVSAPSGIGRLPL LNPIVDHHHHHH
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only