

Recombinant Human Palmitoyl-Protein Thioesterase 1,PPT1 (C-6His)



Catalog No: #AP75563

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75563-1 10ug #AP75563-2 50ug #AP75563-3 500ug #AP75563-4 1mg

Description

Product Name	Recombinant Human Palmitoyl-Protein Thioesterase 1,PPT1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P50897
Uniprot	P50897
GeneID	5538;
Calculated MW	32.3kDa
Target Sequence	DPPAPLPLVIWHGMGDSCCNPLSMGAIKKMVEKKIPGIYVLSLEIGKTLMEDVENSFFLVNSQVTTVCQALAK DPKLQQGYNAMGFSQGGQFLRAVAQRCPSPPMINLISVGGQHQQGVFGLPRCPGESSHICDFIRKTLNAGAYS KVVQERLVQAEYWHDPIKEDVYRNHSIFLADINQERGINESYKKNLMALKKFVMVKFLNDSIVDPVDSEWFGF YRSGQAKETIPLQETSLYTQDRLGLKEMDNAGQLVFLATEGDHLQLSEEWFYAHIIIPFLGVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mM NaCl,10%Glycerol,pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only