

Recombinant Human Probable Serine Carboxypeptidase CPVL,VCP-Like Protein (C-6His)

Catalog No: #AP75564

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75564-1 10ug #AP75564-2 50ug #AP75564-3 500ug #AP75564-4 1mg

Description

Product Name	Recombinant Human Probable Serine Carboxypeptidase CPVL,VCP-Like Protein (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9H3G5
Uniprot	Q9H3G5
GeneID	54504;
Calculated MW	52.9kDa
Target Sequence	<p>LFRSLYRSVSMPPKGDSSGQPLFLTPYIEAGKIQKGRELSLVGPPFGLNMKSYAGFLT VNKTYNSNLFFWFFPA</p> <p>QIQPEDAPVVLWLQGGPGGSSMFGFLFVEHGPYVVTSNMTLRDRDFPWTTL SMLYIDNPVGTGFSFTDDTHG</p> <p>YAVNEDDVARDLYSALIQQFFQIFPEYKNNDFYVTGESYAGKYVPAIAHLIHS LN PVREVKINLNGIAIGDGYSDPE</p> <p>SIIGGYAEFLYQIGLLDEKQKKYFQKQCHECIEHIRKQNWFEAFEILDKLLDGLTSDPSYFQNV TGCSNYYNFL</p> <p>RCTEPEDQLYYVKFSLPEVRQAIHVG NQTFNDGTIVEKYLREDTVQSVK P WLTEIMN NYK VLIYNGQLDIIVAA</p> <p>ALTERSLMGMDWKGSQEYKKA EKKVWKIFKSDSEVAGYIRQV GDFHQVIIRGGGHILPYDQPLRAFDMINRFI</p> <p>YGKGWDPYVGVDHHHHHHH</p>
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only