

Recombinant Mouse Tissue Inhibitor of Metalloproteinases 2, TIMP-2 (C-6His)



Catalog No: #AP75567

Orders: order@signalwayantibody.com

Package Size: #AP75567-1 10ug #AP75567-2 50ug #AP75567-3 500ug #AP75567-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Tissue Inhibitor of Metalloproteinases 2, TIMP-2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P25785
Uniprot	P25785
Calculated MW	22.5kDa
Target Sequence	CSCSPVHPQQAFCNADVIRAKAVSEKEVDSGNDIYGNIPIKRIQYEIKQIKMFKGPKDKDIEFIYTAPSSAVCGVS LDVGGKKEYLIAGKAEGDGKMHITLCDFIVPWDTLSITQKKSLNHRYQMGCECKITRCPMIPCYISSPDECLWM DWWTEKSINGHQAKFFACIKRSDGSCAWYRGAAPPKQEFLDIEDPHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only