

# Recombinant Human TGOLN2,TGN38 Homolog,TGN46(C-6His)



Catalog No: #AP75568

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75568-1 10ug #AP75568-2 50ug #AP75568-3 500ug #AP75568-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TGOLN2,TGN38 Homolog,TGN46(C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O43493
Uniprot	O43493
GeneID	10618;
Calculated MW	38.3kDa
Target Sequence	<p>ATESVKQEEAGVRPSAGNVSTHPSLSQRPGGSTKSHPEPQTPKDSPSKSSAEAQTPEDTPNKSGAEAKTQK  DSSNKSGAEAKTQKGSTSKSGSEAQTTKDSTKSHPELQTPKDSTGKSGAEAQTPEDSPNRSAGAEAKTQKD  SPSKSGSEAQTTKDVPNKSGADGQTPKDGSSKSGAEDQTPKDVPNKSGAEKQTPKDGSNKSGAEEQGPID  GPSKSGAEEQTSKDSPNKVVPEQPSRKDHSPISNPSDNKELPKADTNQLADKGLSPHAFKTESGEETDLIS  PPQEEVKSSEPTEDVEPKEAEDDDTGPEEGSPPKEEKEKMSGASSENREGT LSDSTGSEKDDLYPNGSGN  GSAEVDHHHHHH</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only