

Recombinant Mouse Interleukin-7,IL-7 (C-6His)

Catalog No: #AP75570



Package Size: #AP75570-1 10ug #AP75570-2 50ug #AP75570-3 500ug #AP75570-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Interleukin-7,IL-7 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P10168
Uniprot	P10168
GeneID	16196;
Calculated MW	15.9kDa
Target Sequence	ECHIKDKEGKAYESVLMISIDELDKMTGTDSNCPNNEPNFFRKHVCDTKEAAFLNRAARKLKQFLKMNISEEF NVHLLTVSQGTQTLVNCTSKEEKNVKEQKKNDACFLKRLREIKTCWNKILKGSIVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only