

Recombinant Human GM-CSF,CSF2 (C-6His, Human Cells)

Catalog No: #AP75576

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75576-1 10ug #AP75576-2 50ug #AP75576-3 500ug #AP75576-4 1mg

Description

Product Name	Recombinant Human GM-CSF,CSF2 (C-6His, Human Cells)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P04141
Uniprot	P04141
GeneID	1437;
Calculated MW	15.5kDa
Target Sequence	APARSPSPSTQPWEHVNAIQEARLLNLSRDAAEMNETVEVISEMFDLQEPTCLQTRLELYKQGLRGS�TKL KGPLTMMASHYKQHCPTPETSCATQITFESFKENLKDFLLVIPFDCWEPVQEVDDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only