

Recombinant Human Secretogranin-2, SCG2 (C-6His)

Catalog No: #AP75579



Package Size: #AP75579-1 10ug #AP75579-2 50ug #AP75579-3 500ug #AP75579-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Secretogranin-2, SCG2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P13521
Uniprot	P13521
GeneID	7857;
Calculated MW	68.8kDa
Target Sequence	<p> SFQRNQLLQKEPDLRLLENVQKFPSPEMIRALEYIENLRQQAHKKESSPDYNPYQGVSVPLQQKENGDESHLP ERDSLSEEDWMRIILEALRQAENEPQSAPKENKPYALNSEKNFPMDSDDYETQQWPERKCLKHMQFPMPYE ENSRDNPFKRTNEIVVEEQYTPQSLATLESVFQELGKLTGPNNQKRERMDEEQKLYTDDDDIYKANNIAYEDV VGGEDWNPVEEKIESQTQEEVRDSKENIEKNEQINDEMKRSGQLGIQEEGLRKEKSDQLSDDVSKVIAYLKRL VNAAGSGRLQNGQNGERATRLFEKPLDSQSIYQLIEISRNLIQIPPEDLIEMLKTGEKPNGSVEPERELDLPVDL DDISEADLDHPDLFQNRMLSKSGYPKTPGGAGTEALPDGLSVEDILNLLGMESAAANQKTSYFPNPNQEKVLP RLPYGAGRSRSNQLPKAAWIPHVNRQMAYENLNDKDQELGEYLARMLVKYPEIINSNQVKRVPGQGSSEG DLQEEEQIEQAIKEHLNQGSSQETDKLAPVSKRFPVGGPPKNDTTPNRQYWDEDLLMKVLEYLNQEKAKEGRE HIAKRAMENMVDHHHHHH </p>
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only