

Recombinant Human Leukocyte Ig-Like Receptor A3,LILRA3,ILT6,CD85e (C-6His)



Catalog No: #AP75580

Orders: order@signalwayantibody.com

Package Size: #AP75580-1 10ug #AP75580-2 50ug #AP75580-3 500ug #AP75580-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Leukocyte Ig-Like Receptor A3,LILRA3,ILT6,CD85e (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8N6C8
Uniprot	Q8N6C8
GeneID	11026;
Calculated MW	46.6kDa
Target Sequence	<p>THVQAGPLPKPTLWAEPGSVITQGSPVTLRCQGSLETQEYHLYREKKTALWITRIPQELVKKGQFPILSITWEH AGRYCCYGSHTVGLSESSDPLELVVTGAYSKPTLSALPSPVVTSGGNVTIQCDSQVAFDGFILCKEGEDHP QCLNSHSHARGSSRAIFSVGPVSPSRRWSYRCYGYDSRAPYVWSLPSDLLGLLVPGVSKKPSLSVQPGPVV APGEKLTFCQCGSDAGYDRFVLYKEWGRDFLQRPGRQPQAGLSQANFTLGPVSRSYGGQYTCSGAYNLSSE WSAPSDPLDILITGQIRARPFLSVRPGPTVASGENVTLLCQSQGGMHTFLLTKEGAADSPLRLKSKRQSHKYQ AEFPMSPV TSAHAGTYRCYGLSSNPYLLTHPSDPLELVSGAAETLSPQNKSDSKAGEVDHHHHHH</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only