

Recombinant Human β -1,4-Galactosyltransferase 3,B4GALT3 (C-6His)

Catalog No: #AP75582

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75582-1 10ug #AP75582-2 50ug #AP75582-3 500ug #AP75582-4 1mg

Description

Product Name	Recombinant Human β -1,4-Galactosyltransferase 3,B4GALT3 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O60512
Uniprot	O60512
GeneID	8703;
Calculated MW	41.5kDa
Target Sequence	RLSALFGRDQGPTFDYSHPRDVYSNLSHLPGAPGGPPAPQGLPYCPERSPLLVGPVSVSFSFPVPSLAEIVER NPRVEPGGRYPAGCEPRSRTAIIVPHRAREHHLRLLLYHLHPFLQRQQLAYGIYVIHQAGNGTFNRAKLLNV GVREALRDEEWDCLFLHDVDLLPENDHNLVCDPRGPRHVAVAMNKFYSLYPQYFGGVSALTPDQYLKM NGFPNEYWGWGGEDDDIATRVLGAMKISRPPTSVGHYKMKVHRGDKGNEENPHRFDLLVRTQNSWTQDG MNSLTYQLLARELGPLYTNITADIGTDPRGPRAPSGPRYPGSSQAFRQEMLQRRPPARPGPLSTANHTALR GSHVDHHHHHH
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM TrisHCl, 150mM NaCl, 2mM MgCl ₂ , 10%Glycerol, pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only