

Recombinant Human Interferon ω -1,IFNW1 (C-6His)

Catalog No: #AP75583



Package Size: #AP75583-1 10ug #AP75583-2 50ug #AP75583-3 500ug #AP75583-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Interferon ω -1,IFNW1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P05000
Uniprot	P05000
GeneID	3467;
Calculated MW	21.1kDa
Target Sequence	LGCDLPQNHGLLSRNTLVLLHQMRRISPFLCLKDRRDFRFPQEMVKGSQLQKAHVMSVLHEMLQQIFSLFHTE RSSAAWNMTLLDQLHTGLHQQLQHLETCLLQVVGEGESAGAISSPALTLRRYFQGIRVYLKEKKYSDCAWEV VRMEIMKSLFLSTNMQERLRSKDRDLGSSVDHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB,150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only