

Recombinant Human Follistatin 288,FST (C-6His)

Catalog No: #AP75600



Package Size: #AP75600-1 10ug #AP75600-2 50ug #AP75600-3 500ug #AP75600-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Follistatin 288,FST (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P19883
Uniprot	P19883
GeneID	10468;
Calculated MW	31.5kDa
Target Sequence	GNCWLRQAKNGRCQVLYKTELSKEECCSTGRLSTSWTEEDVNDNTLFKWMIFNGGAPNCIPCKETCENVDC GPGKKCRMNKKNKPRCVCAPDCSNITWKGPCGLDGKTYRNECALLKARCKEQPELEVQYQGRCKKTCRD VFPCGSSTCVVDQTNAYCVTCNRICPEPASSEQYLCGNDGVTYSSACHLRKATCLLGRSIGLAYEGKCIKAK SCEDIQCTGGKKCLWDFKVGRGRCSLDELCPDSKSDEPVCASDNATYASECAMKEAACSSGVLLEVKHSG SCN
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Note: This product is for in vitro research use only