

Recombinant Mouse IL-17A

Catalog No: #AP75607



Package Size: #AP75607-1 10ug #AP75607-2 50ug #AP75607-3 500ug #AP75607-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse IL-17A
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q62386
Uniprot	Q62386
GeneID	16171;
Calculated MW	16.2kDa
Target Sequence	NQGSVAGSCSDRTISSGTQIPQGTLDHIRKYLKAFHRCPPFFIRFQLQSKSVCGGSQDQWVRELVDCFERKE CGTGHGKSFHHQKHLQPASTQTPEAAEGTPSDTSTPAHSQSTQHSTLPSGALSLNKEHTQPWEMTTLPSGY GLEARPEAEANEKQDDRQEQEAPGAGASTPAWVDHHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only