

Recombinant Mouse VEGF-D,PIGF (C-6His)

Catalog No: #AP75611



Package Size: #AP75611-1 10ug #AP75611-2 50ug #AP75611-3 500ug #AP75611-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse VEGF-D,PIGF (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P97946
Uniprot	P97946
GeneID	14205;
Calculated MW	13.1kDa
Target Sequence	FYDTETLKVIDEWQRTQCSPRETCVEVASELGKTTNTFFKPPCVNVFRCGGCCNEEGVMCMNTSTSYISKQ LFEISVPLTSVPELVPVKIANHTGCKCLPTGPRHPYSVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only