

# Recombinant Human Growth Hormone-Releasing Factor, GHRF, Somatoliberin (C-Fc)



Catalog No: #AP75612

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75612-1 10ug #AP75612-2 50ug #AP75612-3 500ug #AP75612-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Growth Hormone-Releasing Factor, GHRF, Somatoliberin (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01286
Uniprot	P01286
GeneID	2691;
Calculated MW	33.4kDa
Target Sequence	<p>PPPPLTLRMRRYADAIFTNSYRKVLGQLSARKLLQDIMSRQQGESNQERGARARLVDDIEGRMDEPKSCDKT  HTCPPCPAPPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQ  YNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCL  VKGFPYSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCSCVMHEALHNHYTQK  SLSLSPGK</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only