

Recombinant Mouse Cystatin F,CST7 (C-6His)

Catalog No: #AP75613



Package Size: #AP75613-1 10ug #AP75613-2 50ug #AP75613-3 500ug #AP75613-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Cystatin F,CST7 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	O89098
Uniprot	O89098
GeneID	13011;
Calculated MW	15.4kDa
Target Sequence	ARPPDFCSKDLISSVKPGFPKTIETNNPGVLKAAARHSVEKFNNDIFLFKESHVSKALVQVVKGLKYMLEVKI GRTTCRKTMMHQLDNCDFQTNPALKRRTLYCYSEVWVWIPWLHSFEVPVLLCQVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only