

Recombinant Human TWEAK Receptor, TWEAK R, TNFRSF12A (C-Fc)



Catalog No: #AP75619

Orders: order@signalwayantibody.com

Package Size: #AP75619-1 10ug #AP75619-2 50ug #AP75619-3 500ug #AP75619-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TWEAK Receptor, TWEAK R, TNFRSF12A (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NP84
Uniprot	Q9NP84
GeneID	51330;
Calculated MW	35.6kDa
Target Sequence	EQAPGTAPCSRGSSWSADLDKCMDCASCRRARPHSDFCLGCAAAPPAPFRLLWVDDIEGRMDEPKSCDKTHT CPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYN STYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVK GFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFSCSVMHEALHNHYTQKSL SLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only