

Recombinant Human Intercellular Adhesion Molecule 1, ICAM-1, CD54 (C-Fc)



Catalog No: #AP75626

Orders: order@signalwayantibody.com

Package Size: #AP75626-1 10ug #AP75626-2 50ug #AP75626-3 500ug #AP75626-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Intercellular Adhesion Molecule 1, ICAM-1, CD54 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P05362
Uniprot	P05362
GeneID	3383;
Calculated MW	76.8kDa
Target Sequence	<p>NAQTSVSPSKVILPRGGSVLVTCTSCDQPKLLGIETPLPKKELLPGNNRKVYELSNVQEDSQPMCYSNCPD</p> <p>GQSTAKFTLVYWTPELAPLPSWQPVGKNLTLRCQVEGGAPRANLTVLLRGEKELKREPAVGEPAEVT</p> <p>TTVLVRRDHHGANFSCRTDLRPQGLELFENTSAPYQLQTFVLPATPPQLVSPRVLEVDTQGTWVCSLDGLF</p> <p>PVSEAQVHLALGDQRLNPTVYGNDSFSAKASVSVTAEDEGTQRLTCAVILGNQSQETLQTVTIYSFPAPNVIL</p> <p>TKPEVSEGTEVTVKEAHPRAKVTLNGVPAQPLGPRAQLLLKATPEDNGRSFSCSATLEVAGQLIHKNTREL</p> <p>RVLYGPRLDERDCPGNWTWPENSQQTPMCQAWGNPLPELKCLKDGTFFPLIGESVTVTRDLEGTYLCRARS</p> <p>TQGEVTRKVTNVLSPRYEVDIEGRMDEPKSCDKTHTCPPCPAPPELLGGPSVFLFPPKPKDTLMISRTPEVT</p> <p>CVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPI</p> <p>EKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFF</p> <p>LYSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only