

Akt1 (Phospho-Thr72) Polyclonal Antibody

Catalog No: #12312

Package Size: #12312-1 50ul #12312-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Akt1 (Phospho-Thr72) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-Akt1 (T72) Polyclonal Antibody detects endogenous levels of Akt1 protein only when phosphorylated at T72.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Akt1 around the phosphorylation site of T72.
Target Name	Akt1
Modification	Phospho
Other Names	AKT1; PKB; RAC; RAC-alpha serine/threonine-protein kinase; Protein kinase B; PKB; Protein kinase B alpha; PKB alpha; Proto-oncogene c-Akt; RAC-PK-alpha
Accession No.	Swiss-Prot: P31749NCBI Gene ID: 207
Uniprot	P31749
GeneID	207;
Target Species	human
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1/500 - 1/2000


Immunohistochemistry: 1/100 - 1/300

ELISA: 1/5000

Not yet tested in other applications

Images

Western Blot analysis of 3T3 cells using Phospho-Akt1 (T72)
Polyclonal Antibody



Note: This product is for in vitro research use only