

Recombinant Human Lymphotoxin β R,LTBR,TNFRSF3,TNFRrp (C-Fc)



Catalog No: #AP75628

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75628-1 10ug #AP75628-2 50ug #AP75628-3 500ug #AP75628-4 1mg

Description

Product Name	Recombinant Human Lymphotoxin β R,LTBR,TNFRSF3,TNFRrp (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P36941
Uniprot	P36941
GeneID	4055;
Calculated MW	48.8kDa
Target Sequence	QAVPPYASENQTCRDQEKEYEPQHRICCSRCPPGTYSVSAKCSRIRDVTCATCAENSYNEHWNLYLTICQLCR PCDPVMGLEEIAPCTSKRKTQCRCQPGMFCAAWALECTHCELLSDCPPGTEAEKDEVGKGNHCVCKAG HFQNTSSPSARCQPHTRCENQGLVEAAPGTAQSDTTCKNPLELPPMSGTMLMVDDIEGRMDEPKSCDKT HTCPPCPAPPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQ YNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCL VKGFYPSDIAVEWESNGQPENNYKTTTPVLDSGDGSFFLYSKLTVDKSRWQQGNVFCSSVMHEALHNHYTQK SLSLSPGK
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only