

# Recombinant Human Cadherin-6,K-Cadherin,CDH6 (C-Fc)



Catalog No: #AP75629

Orders: order@signalwayantibody.com

Package Size: #AP75629-1 10ug #AP75629-2 50ug #AP75629-3 500ug #AP75629-4 1mg

Support: tech@signalwayantibody.com

## Description

Product Name	Recombinant Human Cadherin-6,K-Cadherin,CDH6 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P55285
Uniprot	P55285
GeneID	1004;
Calculated MW	93.5kDa
Target Sequence	TLSTPLSKRTSGFPAKKRALELSGNSKNELNRSKRSMWNQFFLEEYTGSDYQYVVGKLSHQDQRDGLSLK YILSGDGAGDLFIINENTGDIQATKRLDREEKPVYILRAQAINRRTGRPVEPESEFIIKIHDINDNEPIFTKEVYAT VPEMSDVGTFFVQVTATDADDPTYGNSAKVVYSILQGQPYFSVESETGIKTALLNMDRENREYQVVIQAKD MGGQMGGLSGTTTTVNITLTDVNDNPPRFPQSTYQFKTPESPPTPIGRIKASDADVGENAEIEYSITDGEGL DMFDVITDQETQEGIITVKLLDFEKKKVYTLKVEASNPYVEPRFLYLGPFKDSATVRIVVEDVDEPPVFSKLAYI LQIREDAQINTTIGSVTAQDPDAARNPVKYSVDRHTMDRIFNIDSGNGSIFTSKLLDRELLWHNITVIATEINN PKQSSRVPLYIKVLDVNDNAPEFAEFYETFVCEKAKADQLIQLTHAVDKDDPYSGHQFSFLAPEAASGSNFTI QDNKDNTAGILTRKNGYNRHEMSTYLLPVVISDNDYPVQSSTGTVTVRVCACDHHGNMQSCHAEALIHPTGL STGAVDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKF NWWYVDGVEVHNAKTKPREEQYNSTYRVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQ VYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLTVDKSRWQQ GNVFSCSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only