

Recombinant Human NKG2D Ligand 2,NKG2DL2,N2DL2 (C-Fc)

Catalog No: #AP75632

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75632-1 10ug #AP75632-2 50ug #AP75632-3 500ug #AP75632-4 1mg

Description

| | |
|-----------------|---|
| Product Name | Recombinant Human NKG2D Ligand 2,NKG2DL2,N2DL2 (C-Fc) |
| Host Species | Human Cells |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | Q9BZM5 |
| Uniprot | Q9BZM5 |
| GeneID | 80328; |
| Calculated MW | 51.4kDa |
| Target Sequence | GRADPHSLCYDITVIPKFRPGPRWCAVQGGQVDEKTFLLHYDCGNKTVTPVSPLGKKNVTTAWKAQNPVLREV VDILTEQLRDIQLENYTPKEPLTLQARMSCEQKAEGHSSGSWQFSFDGQIFLLFDSEKRMWTTVHPGARKMK EKWENDKVVAMSFHYFSMGDCIGWLEDFLMGMDSTLEPSAGAPLAMSSVDDIEGRMDEPKSCDKTHTCPPC PAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYR VVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYP SDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSP GK |
| Formulation | Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4. |
| Storage | Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months. |

Note: This product is for in vitro research use only