

Recombinant Mouse Cystatin E,CST6 (C-6His)

Catalog No: #AP75633



Package Size: #AP75633-1 10ug #AP75633-2 50ug #AP75633-3 500ug #AP75633-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Cystatin E,CST6 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q9D1B1
Uniprot	Q9D1B1
Calculated MW	14.8kDa
Target Sequence	ELRSRRTGERQNLSPDPRVQKAAQAAVASYNMGSDSLYYFRDTKVIDAKYQLVAGIKYYLTLDIESTECRKT RVSGEHMDLTTCPAAGGQKEKLRNCFELLEVPWKNTTQLLKHDCVQVVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM MES, 150mM NaCl, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only