

Recombinant Human Fas, TNFRSF6, CD95 (C-Fc)

Catalog No: #AP75634



Package Size: #AP75634-1 10ug #AP75634-2 50ug #AP75634-3 500ug #AP75634-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Fas, TNFRSF6, CD95 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P25445
Uniprot	P25445
GeneID	355;
Calculated MW	44.6kDa
Target Sequence	QVTDINSGLELRKTVTTVETQNLLEGLHHDGQFCHKPCPPGERKARDCTVNGDEPDCVPCQEGKEYTDKAH FSSKCRRCRLCDEGHGLEVEINCTRTQNTKCRCKPNFFCNSTVCEHCDCPCTKCEHGIIECTLTSTNKCKEEG SRSNVDDIEGRMDEPKSCDKTHTCPPCPAPPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKF NWWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQ VYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQ GNVFSCSVMHEALHNHYTQKSLSLSPGKHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only