

Recombinant Human Sclerostin,SOST (C-6His)

Catalog No: #AP75636



Package Size: #AP75636-1 10ug #AP75636-2 50ug #AP75636-3 500ug #AP75636-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Sclerostin,SOST (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9BQB4
Uniprot	Q9BQB4
GeneID	50964;
Calculated MW	22.3kDa
Target Sequence	QGWQAFKNDATEIIPELGEYPEPPPELENNKTMNRAENGGRPPHHPFETKDVSEYSCRELHFTRYVTDGPCR SAKPVTELVCSGQC GPARLLPNAIGRGKWWRPSGPDFRCIPDRYRAQRVQLLCPGGEAPRARKVRLVASCK CKRLTRFHNQSELKDFGTEAAPQKGRKPRPRARSAKANQAELENAYHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 250mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only