

Recombinant Human α -Lactalbumin, LALBA (C-6His)

Catalog No: #AP75637



Package Size: #AP75637-1 10ug #AP75637-2 50ug #AP75637-3 500ug #AP75637-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human α -Lactalbumin, LALBA (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P00709
Uniprot	P00709
GeneID	3906;
Calculated MW	15.1kDa
Target Sequence	KQFTKCELSQLLKDIDGYGGIALPELICTMFHTSGYDTQAIVENNESTEYGLFQISNKLWCKSSQVPQSRNICDI SCDKFLDDDDITDDIMCAKKILDIKGIDYWLAHKALCTEKLEQWLCEKLVDDHHHHHH
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM Tris, 150mM NaCl, PH8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only