

# Recombinant Human N-Acetylgalactosaminyltransferase 7, GALNT7 (C-6His)

Catalog No: #AP75638

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75638-1 10ug #AP75638-2 50ug #AP75638-3 500ug #AP75638-4 1mg

## Description

Product Name	Recombinant Human N-Acetylgalactosaminyltransferase 7, GALNT7 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q86SF2
Uniprot	Q86SF2
GeneID	51809;
Calculated MW	73.2kDa
Target Sequence	PRPDDPSPLSRMREDRDVNDPMPNRGGNGLAPGEDRFKPVVPWPHVEGVEVDLESIRRINKAKNEQEHHAG GDSQKDIMQRQYLTFKPQFTFYHDPVLRPGILGNFEPKEPEPPGVVGGPGEKAKPLVLGPEFKQAIQASIKEF GFNMVASDMISLDRSVNDLRQEECKYWHYDENLLTSSVVIVFHNEGWSTLMRTVHSVIKRTPRKYLAEIVLIDD FSNKEHLKEKLDEYIKLWNLVKVFRNERREGLIQARSIGAQAQKAKLGQVLIYLDACEVAVNWYAPLVAPISKD RTICTVPLIDVINGNTYEIIPQGGGDEEDGYARGAWDWSMLWKRVP LTPQEKR LRKTKTEPYRSPAMAGGLFAI EREFFELGLYDPGLQIWGGENFEISYKIWQC CGKLLFVPCSRVGHYRLEGWQGNPPPIYVGSSTLKNYVR VVEVWWDEYKDYFYASRPESQALPYGDISELK KFREDHNCKSFKWFMEEIAYDITSHYPLPPKNVDWGEIRG FETAYCIDSMGKTNGGFVELGPCHRMGGNQLFRINEANQLMQYDQCLTKGADGSKVMITHCNLNEFKEWQY FKNLHRFTHIPSGKCLDRSEVLHQVFISNCDSSKTTQKWEMNNIHSVVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 50mM TrisHCl, 10mM reduced Glutathione, PH8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only