

Recombinant Human PI3K-Interacting Protein 1, PIK3IP1 (C-6His)

Catalog No: #AP75639

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75639-1 10ug #AP75639-2 50ug #AP75639-3 500ug #AP75639-4 1mg

Description

Product Name	Recombinant Human PI3K-Interacting Protein 1, PIK3IP1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96FE7
Uniprot	Q96FE7
GeneID	113791;
Calculated MW	16.7kDa
Target Sequence	SGGCFWDNGHLYREDQTSPAPGLRCLNWLDAQSGLASAPVSGAGNHSYCRNPDEDPRGPWCYVSGEAGV PEKRPCEDLRCPETTSQALPAFTTEIQEASEGPGADEVQVFAPANALPARSEAAAVQPVIGISQRVRMNSKEK KDLGTVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 50mM TrisHCl, 10mM reduced Glutathione, PH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only