

Recombinant Human ER Resident Protein 44,ERp44,TXNDC4 (C-6His)

Catalog No: #AP75640

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75640-1 10ug #AP75640-2 50ug #AP75640-3 500ug #AP75640-4 1mg

Description

Product Name	Recombinant Human ER Resident Protein 44,ERp44,TXNDC4 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9BS26
Uniprot	Q9BS26
GeneID	23071;
Calculated MW	44.2kDa
Target Sequence	EITSLDTENIDEILNNADVALVNFYADWCRFSQMLHPIFEEASDVIKEEFPNENQVVFARVDCDQHSDIAQRYRI SKYPTLKLFRNGMMMKREYRGQRSVKALADYIRQQKSDPIQEIRDLAIEITTLDRSKRNIIGYFEQKSDNYRVF ERVANILHDDCAFLSAFGDVSKPERYSGDNIIYKPPGHSAPDMVYLGAMTNFDVTYNWIQDKCVPLVREITFEN GEELTEEGLPFLILFHMKEDTESLEIFQNEVARQLISEKGTINFLHADCDKFRHPLLHIQKTPADCPVIAIDSRH MYVFGDFKDVLPGLKQFVFDLHSGKLHREFHHGPDPTDTAPGEQAQDVASSPPESSFQKLAPSEYRYTLL RDVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris HCl10%glycerol,PH7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only