

Recombinant Human α -1-Acid Glycoprotein 2,AGP2,ORM2 (C-6His)

Catalog No: #AP75642

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75642-1 10ug #AP75642-2 50ug #AP75642-3 500ug #AP75642-4 1mg

Description

Product Name	Recombinant Human α -1-Acid Glycoprotein 2,AGP2,ORM2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P19652
Uniprot	P19652
GeneID	5005;
Calculated MW	22.6kDa
Target Sequence	<p>QIPLCANLVPVITNATLDRITGKWFYIASAFRNEEYNKSVQEIQATFFYFTPNKTEDTIFLREYQTRQNQCFYN</p> <p>SSYLNQRENGTVSRYEGGREHVAHLLFLRDTKTLMFGSYLDDEKNWGLSFYADKPETTKEQLGEFYEALDC</p> <p>LCIPRSDVMYTDWKKDKCEPLEKQHEKERKQEEGESVDHHHHHH</p>
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only