

# Recombinant Human GDNF Family Receptor α-1,GFRA1 (C-Fc)

Catalog No: #AP75643

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75643-1 10ug #AP75643-2 50ug #AP75643-3 500ug #AP75643-4 1mg

## Description

Product Name	Recombinant Human GDNF Family Receptor α-1,GFRA1 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P56159
Uniprot	P56159
GeneID	2674;
Calculated MW	72.3kDa
Target Sequence	<p>DRLDCVKASDQCLKEQSCSTKYRTLQCVAGKETNFSLASGLEAKDECRSAMEALKQKSLYNCRCKRGMKK  EKNCLRIYWSMYQSLQGNLLEDSPYEPVNSRLSDIFRVVPFISVEHIPKGNCLDAAKACNLDDICKKYRSAYI  TPCTTSVSNVDCNRRKCHKALRQFFDKVPAKHSYGMLFCSCRDIACRERRRQTIVPVCSYEEREKPNCLNLQ  DSCKTNYICRSLADFFTNCPESRSVSSCLKENYADCLLAYSGLIGTVMPNYIDSSSLSVAPWCDCSNSGN  DLEECLKFLNFFKDNTCLKNAIQAFNGSDVTVWQPAFPVQTTTATTTALRVKNKPLGPAGSENIPTHLVLP  CANLQAQKLKSNVSGNTHLCISNGNYEKEGLGASSHITTKVDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPS  VFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQ  DWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESN  GQPENNYKTTTPVLDSGDSFFLYSKLTVDKSRWQQGNVFCSSVMHEALHNHYTQKSLSLSPGK</p>
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only