

# Recombinant Mouse Interleukin-13,IL-13 (Ser26-Phe131,C-6His)

Catalog No: #AP75645

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75645-1 10ug #AP75645-2 50ug #AP75645-3 500ug #AP75645-4 1mg

## Description

Product Name	Recombinant Mouse Interleukin-13,IL-13 (Ser26-Phe131,C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P20109
Uniprot	P20109
GeneID	16163;
Calculated MW	12.7kDa
Target Sequence	SVSLPLTLKELIEELSNITQDQTPLCNGSMVWSVDLAAGGFCVALDSL TNISNCNAIYRTQRILHGLCNRKAPTT VSSLPDTKIEVAHFITKLLSYTKQLFRHGPFVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only