

Recombinant Human Ephrin-A1,EFNA1,LERK1 (C-Fc)

Catalog No: #AP75651



Package Size: #AP75651-1 10ug #AP75651-2 50ug #AP75651-3 500ug #AP75651-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Ephrin-A1,EFNA1,LERK1 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P20827
Uniprot	P20827
GeneID	1942;
Calculated MW	46.5kDa
Target Sequence	DRHTVFWNSSNPKFRNEDYTIHVQLNDYVDIICPHYEDHSVADAAMEQYILYLVEHEEYQLCQPQSKDQVRW QCNRPSTAKHGPEKLSEKQRFPTFTLGKEFKEGHSYYYISKPIHQHEDRCLRLKVTVSGKITHSPQAHVNPQE KRLAADDPEVRVLHSHIAHSVDDIEGRMDEPKSCDKTHTCPPCPAPPELLGGPSVFLFPPKPKDTLMISRTPEVTC VVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIE KTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFL YSKLTVDKSRWQQGNVFCSSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only