

Recombinant Human Tissue-Type Plasminogen Activator, PLAT (C-6His)



Catalog No: #AP75653

Orders: order@signalwayantibody.com

Package Size: #AP75653-1 10ug #AP75653-2 50ug #AP75653-3 500ug #AP75653-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Tissue-Type Plasminogen Activator, PLAT (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P00750
Uniprot	P00750
GeneID	5327;
Calculated MW	61.65kDa
Target Sequence	SYQVICRDEKTKMIYQQHQSWLRPVLRSNRVEYCWNSGRAQCHSVPVKSCSEPRCFNGGTCQQALYFSD FVCQCPEGFAGKCCCEIDTRATCYEDQGISYRGTWSTAESGAECTNWNSSALAQKPYSGRRPDAIRLGLGNH NYCRNPDRDSKPWCYVFKAGKYSSEFCSTPACSEGNSDCYFGNGSAYRGTHSLTESGASCLPWNSMILIGK VYTAQNPSAQALGLGKHNYCRNPDGDAKPWCHVLKNNRRLTWEYCDVPS CSTCGLRQYSQPQFRIGKGLFAD IASHPWQAAIFAKHRRSPGERFLCGGILISSCWILSAAHCFQERFPPHHLTIVILGRTYRVVPGEEEQKFEVEKYI VHKEFDDDTYDNDIALQLKSDSSRCAQESSVVRTVCLPPADLQLPDWTECELSGYGKHEALSPFYSERLKEA HVRLYPSSRCTSQHLLNRTVTDNMLCAGDTRSGGPQANLHDACQGDSSGGPLVCLNDGRMTLVGIISWGLGC GQKDVPGVYTKVTNYLDWIRDNMRPVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM MES, 150mM NaCl, 0.2mM GaCl ₂ , pH5.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only