

Recombinant Human IgG1 Fc(Asp239Glu, Leu241Met)



Catalog No: #AP75665

Orders: order@signalwayantibody.com

Package Size: #AP75665-1 10ug #AP75665-2 50ug #AP75665-3 500ug #AP75665-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IgG1 Fc(Asp239Glu, Leu241Met)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01857
Uniprot	P01857
Calculated MW	25.9kDa
Target Sequence	<p>GSGSDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAK</p> <p>TKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKN</p> <p>QVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNV FSCVMHEALH</p> <p>NHYTQKSLSLSPGK</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only