

Recombinant Human NGAL,Lipocalin-2,LCN2 (C-6His, Human Cells)

Catalog No: #AP75674

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75674-1 10ug #AP75674-2 50ug #AP75674-3 500ug #AP75674-4 1mg

Description

Product Name	Recombinant Human NGAL,Lipocalin-2,LCN2 (C-6His, Human Cells)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P80188
Uniprot	P80188
GeneID	3934;
Calculated MW	21.6kDa
Target Sequence	QDSTSDLIPAPPLSKVPLQQNFQDNQFQGKWWVVLGNAILREDKDPQKMYATIIYELKEDKSYNVTSVLFRK KKCDYWIRTFVPGCQPGEFTLGNIKSYPGLTSLVLRVSTNYNQHAMVFFKKVSNREYFKITLYGRTKELTS ELKENFIRFSKSLGLPENHIVFPVPIDQCIDGVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of PBS, 50% glycerol,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only