

Recombinant Human Glutamine Synthetase, GLUL (C-6His)



Catalog No: #AP75676

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75676-1 10ug #AP75676-2 50ug #AP75676-3 500ug #AP75676-4 1mg

Description

Product Name	Recombinant Human Glutamine Synthetase, GLUL (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P15104
Uniprot	P15104
GeneID	2752;
Calculated MW	43.13kDa
Target Sequence	TTSASSHLNKGIGKQVYMSLPQGEKVQAMYIWDGTGEGLRCKTRTLTLDSEPKCVEELPEWNFDFGSSTLQSEGS NSDMYLVPAAMFRDPFRKDPNKLVLCEVFKYNRRPAETNLRHTCKRIMDMVSNQHPWFGMEQEYTLMGTDG HPFGWPSNGFPGPQGPYYCGVGADRAYGRDIVEAHYRACLYAGVKIAGTNAEVMPAQWFEQIGPCEGISMG DHLWVARFILHRVCEDFGVIATFDPKPIPGNWNGAGCHTNFSTKAMREENGLKYIEEAIEKLSKRHQYHIRAYD PKGGLDNARRLTGFHETSININDFSAGVANRSASIRIPRTVGQEKKGYFEDRRPSANCDPFSVTEALIRTCLLNE TGDEPFQYKNLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 200mM NaCl, 50mM Imidazole, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only