

Recombinant Human PBLD,MAWBP (N-6His)

Catalog No: #AP75677



Package Size: #AP75677-1 10ug #AP75677-2 50ug #AP75677-3 500ug #AP75677-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PBLD,MAWBP (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P30039
Uniprot	P30039
GeneID	64081;
Calculated MW	33.95kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMKLPIFIADAFTARAFRGNPAAVCLLENELDEDMHQKIAREMNLSETAFIRKL HPTDNFAQSSCFGLRWFTPASEVPLCGHATLASAAVLFHKIKNMNSTLTFVTLSGELRARRAEDGIVLDLPLYP AHPQDFHEVEDLIKTAIGNTLVQDICYSPDTQKLLVRLSDVYNRSFLENLKVNTENLLQVENTGKVKGLILTLKG EPGGQTQAFDFYSRYFAPWVGVAEDPVTGSAHAVLSSYWSQHLGKKEMHAFQCShrGGELGISLRPDGRV DIRGGAAVVLEGLTA
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 100mM NaCl, 1mM DTT, 30% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only