

Recombinant Human Small Ubiquitin-Related Modifier 3,SUMO3,SMT3A (N-6His)



Catalog No: #AP75678

Orders: order@signalwayantibody.com

Package Size: #AP75678-1 10ug #AP75678-2 50ug #AP75678-3 500ug #AP75678-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Small Ubiquitin-Related Modifier 3,SUMO3,SMT3A (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P55854
Uniprot	P55854
GeneID	6612;
Calculated MW	13.8kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMSEEKPKKEGVKTENDHINLKVAGQDGSVVQFKIKRHTPLSKLMKAYCERQG LSMRQIRFRFDGQPINETDTPAQLEMEDEDTIDVFQQQTGGVPESLAGHSF
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only