

Recombinant Human 3-Hydroxybutyrate Dehydrogenase Type 2, BDH2, DHRS6 (N-6His)

Catalog No: #AP75679

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75679-1 10ug #AP75679-2 50ug #AP75679-3 500ug #AP75679-4 1mg

Description

Product Name	Recombinant Human 3-Hydroxybutyrate Dehydrogenase Type 2, BDH2, DHRS6 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9BUT1
Uniprot	Q9BUT1
GeneID	56898;
Calculated MW	28.9kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMGRLDGKVIILTAAQIGIGQAAALAFAREGAKVIATDINESKLQELEKYPGIQT RVLDVTKKKQIDQFASEVERLDVLFNVAGFVHHGTVLDCSEKDWDFSMNLNVRSMYLMIKAFPKMLAQKSG NIINMSSVASSVKGVVNRCVYSTTKAAVIGLTKSVAADFIQQGIRCNCVCPGTVDTPSLQERIQARGNPEEARN DFLKRQKTGRFATAEEIAMLCVYLASDESAYVTGNPVIIDGGWSL
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only