

Recombinant Human Peroxiredoxin-6, PRDX6 (N-6His)

Catalog No: #AP75680

Orders: order@signalwayantibody.com

Package Size: #AP75680-1 10ug #AP75680-2 50ug #AP75680-3 500ug #AP75680-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Peroxiredoxin-6, PRDX6 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P30041
Uniprot	P30041
GeneID	9588;
Calculated MW	27.2kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMPGGLLLGDVAPNFEANTTVGRIRFHDFLGDSWGILFSHPDRFTPVCTTELG RAAKLAPEFAKRNVKLIALSIDSVEDHLAWSKDINAYNCEEPTEKLPFPIIDDRNRELAILLGMLDPAEKDEKGM PVTARVVFVFGPDKKLSILYPATTGRNFDEILRVVISLQLTAEKRVATPVDWKDGDSVMVLPTIPEEEAKKLF PKGVFTKELPSGKKYLRYTPQP
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only