

Recombinant Human Phosphatidylinositol Transfer Protein α Isoform,PITPNA (N-6His)



Catalog No: #AP75683

Orders: order@signalwayantibody.com

Package Size: #AP75683-1 10ug #AP75683-2 50ug #AP75683-3 500ug #AP75683-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Phosphatidylinositol Transfer Protein α Isoform,PITPNA (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q00169
Uniprot	Q00169
GeneID	5306;
Calculated MW	34kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMVLKKEYRVILPVSVDEYQVGQLYSVAEASKNETGGGEGVEVLVNEPYEKD GEKGQYTHKIYHLQSKVPTFVRMLAPEGALNIHEKAWNAYPYCRTVITNEYMKEDFLIKIETWHKPDLTQENV HKLEPEAWKHVEAVYIDIADRSQVLSKDYKAEEDPAKFKSIKTGRGPLGPNWKQELVNQKDCPYMCAYKLVT VKFKWWGLQNKVENFIHKQERRLFTNFHRQLFCWLDKWVDLTMDDIRRMEEETKRQLDEMQRKDPVKGMT ADD
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM TrisHCl, 1mM EDTA, 1mM DTT, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only