

# Recombinant Human Endothelial Differentiation-Related Factor 1,EDF1,MBF1 (C-6His)

Catalog No: #AP75684

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75684-1 10ug #AP75684-2 50ug #AP75684-3 500ug #AP75684-4 1mg

## Description

Product Name	Recombinant Human Endothelial Differentiation-Related Factor 1,EDF1,MBF1 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O60869
Uniprot	O60869
GeneID	8721;
Calculated MW	17.4kDa
Target Sequence	AESDWDTVTVLRKKGPTAAQAKSKQAILAAQRRGEDVETSKKWAAGQNKQHSITKNTAKLDRETEELHHDRV TLEVGKVIQQGRQSKGLTQKDLATKINEKPQVIADYESGRAIPNNQVLGKIERAIGLKLGRGKDIGKPIEKGPRAK LEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only