

Recombinant Human Chloride Intracellular Channel Protein 4, CLIC4 (N-6His)



Catalog No: #AP75686

Orders: order@signalwayantibody.com

Package Size: #AP75686-1 10ug #AP75686-2 50ug #AP75686-3 500ug #AP75686-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Chloride Intracellular Channel Protein 4, CLIC4 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9Y696
Uniprot	Q9Y696
GeneID	25932;
Calculated MW	30.9kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMALSMPLNGLKEEDKEPLIELFVKAGSDGESIGNCPFSQRLFMLWLKGVVF SVTTVDLKRKPADLQNLAPGTHPPFITFNSEVKTDVNIKEEFLEEVLCPKYLKLSPKHPESNTAGMDIFAKFSA YIKNSRPEANEALERGLLKTQLDEYLN SPLPDEIDENSMEDIKFSTRKFLDGNEMTLADCNLLPKLHIVKVA KKYRNFDPKEMTGIWRYLTNAYS RDEFTNTCPSDKEVEIAYS DVAKRLTK
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 1mM DTT, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only